

Identification of the product

Product Name	Gastroenteritis / parasite panel I : Giardia lamblia / Entamoeba histolytica / Cryptosporidium parvum real-time PCR kit		
Catalog Number	DIA-GE/Parasite I-050	Production Date	2015-03
Lot Number	08005C	Expiry Date	2017-03

Analysis procedure

Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	QuantiFast Multiplex PCR + R Kit (Qiagen)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

QC Results by Real Time PCR

Testing	Specifications	Results		
Positive Control (Ct/Cp)	30.0 ± 3.3	FAM : 28.3	YD : 30.1	TR : 29.7
Negative Control (Ct/Cp)	Negative	Negative	Negative	Negative

Approval

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: F. Noye RA officer

Date: 27 May 2015

