

Identification of the product

Product Name	Mycobacterium tuberculosis real-time PCR kit		
Catalog Number	DIA-MT-050	Production Date	2015-07
Lot Number	08005G	Expiry Date	2017-07

Analysis procedure

Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	QuantiFast Multiplex PCR-R kit (Qiagen)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

QC Results by Real Time PCR

Testing	Specifications	Results
Positive Control (Ct/Cp)	30.0 ± 3.3	FAM: 29.0
Negative Control (Ct/Cp)	Negative	Negative

Approval

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: R. ZALETELS

Date: 06 Oct 2015

QA Manager

