

**Identification of the product**

Product Name	R-DiaFlu		
Catalog Number	DDGR-30-L100	Production Date	2014-08
Lot Number	03004H	Expiry Date	2016-08

**Analysis procedure**

Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	TaqMan Fast Virus 1-Step Master Mix (Applied Biosystems)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

**QC Results by Real Time PCR**

Testing	Specifications	Results
Positive Control (Ct/Cp)	FAM : 30.0 ± 3.3	FAM : 28.4
Positive Control (Ct/Cp)	YD : 30.0 ± 3.3	YD : 30.0
Negative Control (Ct/Cp)	Negative	Negative

**Approval**

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: R. ZALTELI  Date: 24 Feb 2016  
*QA/QC Manager*