

Identification of the product

Product Name	R-DiaFlu		
Catalog Number	DDGR-30-L100	Production Date	10-2015
Lot Number	15315J	Expiry Date	04-2017

Analysis procedure

Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	TaqMan Fast Virus 1-Step (Applied Biosystems)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

QC Results by Real Time PCR

Testing	Specifications	Results
Positive Control (Ct/Cp) FluA	30.0 ± 3.3	FAM: 27.4
Positive Control (Ct/Cp) FluB	30.0 ± 3.3	YD: 28.3
Negative Control (Ct/Cp)	Negative	Negative

Approval

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: R. ZALETELS  Date: 12 NOV 2015
QA/QC Manager