

**Identification of the product**

Product Name	R-DiaRes		
Catalog Number	DDGR-46-L100	Production Date	03-2015
Lot Number	15115C	Expiry Date	03-2017

**Analysis procedure**


Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	TaqMan Fast Virus 1-Step (Applied Biosystems)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

**QC Results by Real Time PCR**

Testing	Specifications	Results
Positive Control (Ct/Cp)	30.0 ± 3.3	FAM: 27.7
Positive Control (Ct/Cp)	30.0 ± 3.3	YD: 28.7
Negative Control (Ct/Cp)	Negative	Negative

**Approval**

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: A. ZALTECT  Date: 03 Feb 2016  
*PA (QC Manager)*