

Identification of the product

Product Name	R-DiaRhino		
Catalog Number	DDGR-60-L100	Production Date	10-2015
Lot Number	15305J	Expiry Date	04-2017

Analysis procedure


Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	TaqMan Fast Virus 1-Step (Applied Biosystems)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

QC Results by Real Time PCR

Testing	Specifications	Results
Positive Control (Ct/Cp): HRV	30.0 ± 3.3	FAM: 27.2
Negative Control (Ct/Cp)	Negative	Negative

Approval

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: R. PALETTELS  Date: 19 NOV 2015
QA/QC Manager