

**Identification of the product**

Product Name	S-DiaCT		
Catalog Number	DDGS-10-L100	Production Date	2015-07
Lot Number	12005G	Expiry Date	2017-07

**Analysis procedure**

Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	LightCycler 480 Probes Master (Roche)
Real Time PCR equipment (Supplier)	LightCycler 480 II (Roche)

**QC Results by Real Time PCR**

Testing	Specifications	Results
Positive Control (Ct/Cp)	30.0 ± 3.3	FAM : 30.4
Negative Control (Ct/Cp)	Negative	Negative

**Approval**

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

Approved by: F. Naye QA office

Date: 15 sept 2015

