

Identification of the product

Product Name	S-DiaMGTV		
Catalog Number	DDGS-60-L100	Production Date	2015-06
Lot Number	08005F	Expiry Date	2017-06

Analysis procedure

Volume PCR samples (µL)	5
Volume PCR mix (µL)	25
Master Mix reference (Supplier)	Optima Du Master Mix x5 (Diagenode)
Real Time PCR equipment (Supplier)	LightCycler® 480 II (Roche)

QC Results by Real Time PCR

Testing	Specifications	Results
Positive Control (Ct/Cp) FAM :	30.0 ± 3.3	31.5
Positive Control (Ct/Cp) YD :	30.0 ± 3.3	30.8
Negative Control (Ct/Cp)	Negative	Negative

Approval

This lot has been reviewed by Quality Assurance in compliance of Diagenode's Quality System.

 Approved by: A ZACETELS


 Date: 01 sep 2015

QA/QC Manager